## **PrimeGene** a biotechne brand

## Recombinant Rat Eotaxin/CCL11 (rRtEotaxin/CCL11)

## **PrimeGene Technical Data Sheet**

Catalog Number:	244-11
Source:	Escherichia coli.
Molecular Weight:	Approximately 8.4 kDa, a single non-glycosylated polypeptide chain containing 74 amino acids.
Quantity:	5µg/20µg/1000µg
AA Sequence:	HPGSIPTSCC FTMTSKKIPN TLLKSYKRIT NNRCTLKAIV FKTKLGKEIC ADPKKKWVQD
	ATKHLDQKLQ TPKP
Purity:	> 96 % by SDS-PAGE and HPLC analyses.
<b>Biological Activity:</b>	Fully biologically active when compared to standard. The biological activity determined by a
	chemotaxis bioassay using purified eosinophils is in a concentration range of 0.1-1.0 µg/ml.
Physical Appearance:	Sterile Filtered White lyophilized (freeze-dried) powder.
Formulation:	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 7.4, 150 mM NaCl.
Endotoxin:	Less than 1 EU/µg of rRtEotaxin/CCL11 as determined by LAL method.
<b>Reconstitution:</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the
	bottom. Reconstitute in sterile distilled water or aqueous buffer containing 0.1 % BSA to a
	concentration of 0.1-1.0 mg/mL. Stock solutions should be apportioned into working aliquots and
,	stored at $\leq$ -20 °C. Further dilutions should be made in appropriate buffered solutions.
Shipping:	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature
	recommended below.
Stability & Storage:	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
	• 12 months from date of receipt, -20 to -70 °C as supplied.
	• 1 month, 2 to 8 °C under sterile conditions after reconstitution.
	• 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Usage:	This material is offered by Shanghai PrimeGene Bio-Tech for research, laboratory or further
	evaluation purposes. NOT FOR HUMAN USE.

## Rat Eotaxin/CCL11

CCL11 encoded by the gene CCL11 is belonging to the CC chemokine family. It is a potent eosinophil chemoattractant that was originally purified from bronchoalveolar lavage fluid of guinea pigs sensitized by aerosol challenge with ovalbumin. CCL11 is a strong and specific eosinophil chemoattractant in vitro. It can directly chemotactic for eosinophils, but not for monocytes or neutrophils. CCR3 has been identified to be a specific CCL11 receptor. CCR3 has also been shown to serve as a cofactor for a restricted subset of primary HIV viruses and binding of CCL11 to CCR3 inhibited infection by the HIV isolates.

Website: <u>www.primegene.com</u> Fax: +86 21 61077348